

## PREPARATION OF REACTIVE SOLUTIONS

### PURPOSE AND RATIONALE

To speak to establish a general procedure for preparation reactive solutions.

### MATERIAL

Magnetic Stirrers.  
 Analytical Balance.  
 Top-loading balance.  
 Desiccator.  
 Funnels  
 Drying stove.  
 Filter paper, fast.  
 Polyethylene bottles (of 1 liter) with screw cap.  
 Volumetric flasks of different volumes.  
 Filter paper.  
 Weighing bottle.  
 Pipettes of different volumes.  
 Plates heaters.  
 Rods.  
 Beakers of various volumes.

### REAGENTS

(In addition to the practice specific)  
 Concentrated hydrochloric acid *AG*.  
 Distilled water.  
 Concentrated ammonia *AG*.

### METHODOLOGY

The appropriate methodology varies depending on the characteristics of each reactive solution to prepare.

### CALCULATIONS

Based on the criterion of having data on the concentration of the solution that we want to prepare in grams / liter. To switch from normal to grams / liter applying the expression:

$$c = N \cdot p_e$$

where "c" is the concentration in grams/liter, "N" is the normality and "pe" is the equivalent weight. The amount to be weighed to prepare a volume "v" of concentration "c" is:

$$m = c \cdot v$$

If the amount is too small, we can weigh a higher amount for subsequent dissolution. To prepare a diluted solution from another more concentrated, the volume of take concentrated solution is:

$$V_c = \frac{V_d \cdot C_d}{C_c}$$

where "Vc" is the amount to be taken in concentrated solution, "Vd" the volume of diluted solution to be prepared, "Cd" is the concentration (in grams/liter, molarity or normality) and "Cc" is the concentration of concentrated solution (in grams/liter, molarity or normality).

## OBSERVATIONS

Always follow the precautions of use inherent in each product.

---

### Questionnaire 3.1.- Preparation of reactive solutions

1.- Define the following concepts:

- a) reactive AG
- b) titrated solution
- c) class A volumetric material
- d) class B volumetric material
- e) volume content
- f) volume for spill
- g) pure reagent
- h) extemporaneous solution

2.- Precautions for preparation and preservation of reagents of this characteristics:

- a) highly corrosive solutions of reagents
- b) reducing reagent solutions
- c) solutions of reactive oxidants
- d) very unstable solutions